K-6530-RP CALIBRATION SOLUTIONS (for use with Myron L TechPro II™ TPH1 meter)

COMPONENTS:

1 x R-0833-C	DI Water, 2 oz, Dispenser Cap
--------------	-------------------------------

1 x R-0834-C pH Soaker Solution, 2 oz, Dispenser Tip

1 x R-0868-25C-C Conductivity Solution 2500 μ S (1273 ppm NaCl), 2 oz,

Dispenser Tip

1 x R-1099-04-C Buffer Solution pH 4.0, 2 oz, Dispenser Tip 1 x R-1099-07-C Buffer Solution pH 7.0, 2 oz, Dispenser Tip 1 x R-1099-10-C Buffer Solution pH 10.0, 2 oz, Dispenser Tip

CAREFULLY READ AND FOLLOW PRECAUTIONS ON REAGENT LABELS.
KEEP REAGENTS AWAY FROM CHILDREN.

TO ORDER REPLACEMENT REAGENTS CALL TOLL-FREE 800-TEST KIT (837-8548).

For specific procedures, refer to Operation Manual included with meter.

