

**K-6560-RP  
CALIBRATION SOLUTIONS  
(for use with Myron L EP-10 meter)**

**COMPONENTS:**

- 2 x R-0833-C      DI Water, 2 oz, Dispenser Cap
- 2 x R-0868-25C-C      Conductivity Solution 2500  $\mu$ S (1273 ppm NaCl), 2 oz,  
Dispenser Tip

**CAREFULLY READ AND FOLLOW PRECAUTIONS ON REAGENT LABELS.  
KEEP REAGENTS AWAY FROM CHILDREN.**

**TO ORDER REPLACEMENT REAGENTS CALL TOLL-FREE  
800-TEST KIT (837-8548).**

**For specific procedures, refer to Operation Manual included with meter.**

**Instr. #5848**



31 Loveton Circle, Sparks, MD 21152 USA  
800-TEST KIT (837-8548) • 410-472-4340

111717